

# BIS SHAKERS™

## The Next Generation of Shakers

### Microplate Shaker

#### FEATURES

- Microplate Footprint
- RS 232 Communication Port
- Automation Friendly
- Compatible with microplates, optimized for 384 well plates
- 24 VDC operation (integrated units)
- GUI for stand-alone operation with heaters and other external devices
- Full integration with liquid handling devices

#### APPLICATIONS

- Dissolution Assays
- Plasmid Purification
- Magnetic Bead Assays



**Microplate Shaker**



**Test Tube Shaker**



### Test Tube Shaker

#### FEATURES

- Compact Footprint
- RS 232 Communication Port
- Automation Friendly
- Compatible with Test tubes and Vials
- 24 VDC operation (integrated units)
- GUI for stand-alone operation with heaters and other external devices
- Full integration with liquid handling devices

#### APPLICATIONS

- Dissolution Assays
- Plasmid Purification
- Magnetic Bead Assays




**BioIntegrated**  
SOLUTIONS

Empowering BioTechnology™

## RS 232 Communications Protocol with "Auto Detect"

**BIS Microplate Orbital Shaker**



**Set Parameters**

5 Seconds  
Set Acceleration  
500 RPM  
Set Velocity

☒ CW Rotation  
☐ CCW Rotation

Set Direction

☒ Find Home  
☒ Index After Stopping

☒ Enable Status Text  
☒ Auto-Send

**Motion Commands**

Go  
Stop  
Find Home

**Get Product Info**

Read Product Type  
Get Serial #  
Get Software Version #

**Get Values**

Acceleration  
Velocity  
Current Velocity  
Measured Speed  
Get Text Status

**Diagnostic Tests**

Random Cycling  
Shaft Sensor  
Set Top RPM

**Alternating Speeds**

**Speed #1**

1 Acceleration  
500 RPM  
10 Delay Secs

☒ CW Rotation  
☐ CCW Rotation

**Speed #2**

1 Acceleration  
3500 RPM  
1 Delay Secs

☒ CW Rotation  
☐ CCW Rotation

Number of Cycles 1  
Begin Cycling

**Your RS232 Port for Shakers**

☐ Port 1 disabled/non-existent  
☐ Port 2 disabled/non-existent  
☒ Comm Port 3  
☐ Port 4 disabled/non-existent  
☐ Port 5 disabled/non-existent  
☐ Comm Port 6  
☐ Port 7 disabled/non-existent  
☐ Port 8 disabled/non-existent

Auto-Find Shakers

Software Version 1.0

Instrument Configuration	Bis Microplate Shaker™	Bis Test Tube Shaker™
Dimensions W x D x H	5.500" x 3.500" x 1.500" (140.0 mm x 90.0 mm x 40.0 mm)	5.500" x 8.000" x 5.500" (140.0 mm x 203.0 mm x 140.0 mm)
Power Supply	24 VDC, 180 mA. for integrated units Stand-Alone Units 100 – 240 VAC 50 – 60 Hz	24 VDC, 180 mA. Stand-Alone Units 100 – 240 VAC 50 – 60 Hz
Orbital Motion	2 mm diameter, circular shape, constant everywhere on well plate	4 mm diameter, circular shape, constant everywhere on well plate
Material	Aluminum, nickel plated, stainless steel hardware, white polyethylene clips	Aluminum, nickel plated, stainless steel hardware, white polyethylene clips
Weight	2 lbs	5 lbs
Regulatory Standard	Compliance with CE/CSA	Compliance with CE/CSA
Pricing Information	OEM Packages Available (847) 712-3085; info@biointsol.com	OEM Packages Available (847) 712-3085; info@biointsol.com

**Bio Integrated Solutions Inc**  
**1299 North Arbor Lane, Palatine IL 60067**  
(847) 712-3085 ; (815) 642-0975 (FAX)  
[www.biointsol.com](http://www.biointsol.com)